

A tutorial on spike-field methods made for Vanderbilt neuroscientists

Prof. Dr. Martin Vinck

01.06.2023

Learning goals




- 1) What are the main kinds of spike-field analyses?
- 2) What are the strengths and limitations of different methods?
- 3) How do we practically implement these analyses?

Recommended background literature

Review Article | [Published: 25 June 2018](#)

Investigating large-scale brain dynamics using field potential recordings: analysis and interpretation


[Bijan Pesaran](#) , [Martin Vinck](#), [Gaute T. Einevoll](#), [Anton Sirota](#), [Pascal Fries](#), [Markus Siegel](#), [Wilson Truccolo](#), [Charles E. Schroeder](#) & [Ramesh Srinivasan](#)

[Nature Neuroscience](#) **21**, 903–919 (2018) | [Cite this article](#)

17k Accesses | **180** Citations | **58** Altmetric | [Metrics](#)

Motivation

- Widely used method to study mechanisms and function of neural synchronization
- Highly sensitive method to quantify rhythmicity of different cell types to network activity
- LFP allows to immediately assign a “phase” for every spike relative to network activity
 - LFP picks up activity of about 10^5 neurons in vicinity electrode




Live Drawing



Pipeline




- 1) Measure phase
- 2) Quantify phase locking
- 3) Statistics on phase locking for individual neurons
- 4) Statistics for groups of neurons



Pipeline




- 1) Measure phase
- 2) Quantify phase locking
- 3) Statistics on phase locking for individual neurons
- 4) Statistics for groups of neurons



Methods to estimate phase

- 1) **Short-term Fourier** with fixed window centered on each spike - slow for many neurons
- 2) **Wavelets** (fixed numbers of cycles per frequency) (see Vinck et al., 2013, Neuron)
- 3) **Filtering** -> Hilbert or other forms of phase estimation (see Onorato, ..., Vinck, 2020, Neuron)
- 4) **Other decomposition techniques**, e.g. ICA, empirical mode decomposition



Live Drawing




Live Drawing



Spiketriggeredspectrum with fft



- **Short term Fourier** of LFP centered around each individual spike
 - **Speed trade-off:** Slow if many units and spikes, fast if few and sparse neurons
- **Disadvantage:** The chosen window length is likely only optimal for a narrow frequency range



Live Drawing




Spiketriggeredspectrum with "convol"



- **Wavelet decomposition** with fixed number of cycles per frequency - computed for all t
 - Fast method when data has many units and spikes
- **How to choose number of cycles?**
 - Using too many cycles can underestimate phase locking if oscillation is not narrow-band
 - Using too few cycles brings risk of leakage


Centered on each spike



Pipeline



- 1) Measure phase
- 2) Quantify phase locking
- 3) Statistics on phase locking for individual neurons
- 4) Statistics for groups of neurons




Methods to assess spike-field locking

- **Frequency resolved**

- Spike-field locking
 - Quantify consistency of locking based on set of spike-field phases
- Spike-field coherence
 - Linear correlation between spike trains and LFPs in the frequency domain

- **Time resolved**

- Spike-triggered average



Spike-triggered average (STA)

- **STA:**
 - Take each spike and compute the average LFP around the spike.

Spike-triggered average (STA)

- **STA:**
 - Take each spike and compute the average LFP around the spike.


Spike-triggered average (STA)

- **STA:**

- Take each spike and compute the average LFP around the spike.

- **Properties:**

- STA is not normalized for power - dominated by frequencies with more LFP power
- The STA gives an indication of rhythmicity (side-lobes)
 - Rhythmic side lobes are a property of the LFP signal, not of the spike train




Pairwise phase consistency

- **Phase-locking value - resultant vector length**
 - Biased by spike counts
 - Biased by spiking history effects (refractory period, bursting)

Pairwise phase consistency



- **PPC1: Pairwise phase consistency** (Vinck et al., 2012, J Comp Neuro)
 - Compares pairs of spike phases across trials - like “spin model”
 - Avoids history effect
 - Unbiased by spike count



Live drawing



Spike-field coherence

- **Standard coherence expression** see (Bastos' tutorial)
 - Explained variance for linear prediction at each individual frequency


Spike-field coherence

- **Standard coherence expression** see (Bastos' tutorial)
 - Explained variance for linear prediction at each individual frequency
- **Advantage:**
 - Fast to compute
 - Can deal better with edges of trial periods than spike-centered methods
- **Disadvantages:**
 - Spike-field coherence is a strictly linear measure
 - Not adequate to measure locking to oscillator in broad frequency band
 - Strongly positively biased by firing rates
 - Subsampling approaches are computationally expensive and lose information

Pipeline



- 1) Measure phase
- 2) Quantify phase locking
- 3) Statistics on phase locking for individual neurons
- 4) Statistics for groups of neurons



Single-site statistics

- **Rayleigh test**
 - Circular statistics test for significance individual neurons
 - Does not adequately deal with history effects in spike trains

Single-site statistics

- **Permutation test**

- Preferred statistical testing, permuting spikes and fields across trials
- Can be combined with multiple comparisons

Descriptive group statistics



- Reliability of measurements depends on the number of spikes per neuron

Descriptive group statistics: solutions

- **Unweighted PPC**
 - Can be strongly driven by variance in neurons with few spikes
- **Threshold >50 spikes**
 - Works well in practice, but the threshold is arbitrary
- **Weighted PPC** (Vinck et al., 2013, Spyropoulos et al., 2024, Neuron)
 - Assumes identical distribution of neurons independent of spike count
- **Questionable procedures:**
 - Select neurons with positive PPC values (recreates bias)
 - Select only significantly locked neurons (recreates bias)

Software



- **SPIKE toolbox for Matlab:**

- Main developer *Martin Vinck*

- **Tutorial using Matlab:**

- Simple example of synthetic dataset of spike trains locked to synthetic modulation signal
- Signal is generated as a damped harmonic oscillator driven by noise (Spyropoulos et al.)

- **Optimized Python code** performing similar analyses will be released within ~1 year

Spike structure



- ***spike.time***: {1x#Neurons} cell array

- Each cell contains a 1x#spikes vector: spike times relative to some event

- `spike.time{1} = [0.4360 0.2360 0.7720 0.2520 0.9920 0.1920 0.8000]`

- ***spike.trial***: {1x#Neurons} cell array

- Each cell contains a 1x#spikes vector: the trial in which the spike occurred

- `spike.trial{1} = [3 8 13 16 20 24 24 28 28 28]`

- ***spike.trialtime***: [nTrials x 2] array

- Each row contains the beginning and end of each trial relative to the event

- `spike.trialtime(2,:) = [0 0.9960]`

Spike structure in Matlab

- Working with spike times is much more memory efficient than binary (continuous) format
- For some analyses it is necessary to convert to binary format
 - e.g. Spike-field coherence
- Can convert spike structure to continuous data structure in FieldTrip, and back
- **Important:**
 - Match number of trials and time axis in trial definition of spikes and LFPs

```
>> spike = ft_checkdata(dataspike, 'datatype', 'spike', 'feedback', 'yes');  
the input is raw data with 50 channels and 100 trials  
converting raw data into spike data
```

Spiketriggeredspectrum with fft

- Example using SPIKE toolbox:
 - `cfg.timwin = [-0.1 0.1]`
 - `cfg.taper = 'hann'`
 - `sts = ft_spiketriggeredspectrum(cfg,data,spike)`
- Output added:
 - `sts.fourierspctrm = {1x#Neurons} cell array`
- Each cell contains a
 - `#spikes x #lfpchans x #freq array`

Spiketriggeredspectrum with convol

- Example using SPIKE toolbox with 5 cycles for each frequency and a Hann taper:
 - `cfg.foi = [10 : 1 : 100]`
 - `cfg.taper = 'hann'`
 - `cfg.t_ftimwin = 5 . / cfg.foi;`
 - `sts = ft_spiketriggeredspectrum_convol(cfg,data,spike)`
- Output added:
 - `sts.fourierspctrm = {1x#Neurons} cell array`
- Each cell contains a
 - `#spikes x #lfpchans x #freq array`

Spike-triggered average (STA)

- Example using SPIKE toolbox:
 - `cfg.timwin = [-0.1 0.1]`
 - `sta = ft_spike-triggered-average(cfg,data)`
 - Here data contains both lfp and spiking data in continuous format
- Output added:
 - `sta.avg = #Neurons x #time_points` array
 - `sta.time = 1x#time_points` array